

ESPEC CORP.  
Apr. 2002 revision 2


## **Outline of ERC-100SII**

# Outline of ERC-100SII

## 1. About ERC-100SII

- ERC-100SII is a PC application software for remote monitoring and operating of instrumentation from Windows PC
- ERC-100SII realizes such function by means of PC that the test condition setting, operation status monitoring and operation data recording.
- One PC do the remote monitoring and operation of one instrumentation (1:1).
- RS-232C is used to connect the communication port of PC with the communication port of the back of the chamber.

## 2. System Construction



## 3. Required Equipment

Optional communication board(RS-232C)  
RS-232C cable  
Windows PC